

Population Analysis Report #593-193

Study #3

Report date: 19 July 2019

Prepared For: 3B Medical

799 Overlook Dr

Winter Haven FL 33884-1671

Sample Identification #: Mesa Labs Bacillus pumilus Spore Disc Biological Indicators,

DPSSC/3

Test Article: CPAP Accessories – Lumin CPAP Cleaner

Date Testing Performed: 16-Jul-2019
Date Testing Completed: 18-Jul-2019

Test Procedure: Per Mesa Labs' population assay procedures

Procedure: A population assay was performed on the Mesa Labs Spore Discs following Mesa Labs' procedures. Plates were incubated at 30-35°C and results were recorded after 48 hours. Pictures were taken of the BIs upon receipt at Mesa Labs (Attachment 1) and of the plates after incubation (Attachment 2).

Population Results:

Table I: Exposed Sample Results

Sample ID	Description	Duration	Verified Population	
G1	CPAP Mask	5 Min	0	
G2	CPAP Mask	5 Min	0	
H1	CPAP Hose	5 Min	0	
H2	CPAP Hose	5 Min	0	
I1	CPAP Water Chamber	5 Min	2.5000 x 10°*	
I 2	CPAP Water Chamber	5 Min	0	

^{*}Estimated population only due to less than 6 CFUs.

Table II: Unexposed (Positive Control) Results

Sample ID	Verified Population	
CN3	2.1300×10^{3}	

Mesa	Labor	atories,	Inc.

3B Medical

Population Assays on Exposed an Unexposed DPSSC/3 Biological Indicators

Prepared by:

Samantha Bodnar, Contract Studies Lab Analyst Mesa Labs, Bozeman Manufacturing Facility Date: ZSul 2019

Approved by:

Jared Turner, Contract Studies Lab Coordinator Mesa Labs, Bozeman Manufacturing Facility Date: 25JUL 2019


Approved by

Nicole DeWald, Quality Assurance Manager Mesa Labs, Bozeman Manufacturing Facility Date: 25542019

Population Assays on Exposed an Unexposed DPSSC/3 Biological Indicators

Attachment 1


Sample Receipt Pictures



3B Medical

Population Assays on Exposed an Unexposed DPSSC/3 Biological Indicators


Attachment 1 (cont.)



Attachment 1 (cont.)



Attachment 1 (cont.)




3B Medical

Population Assays on Exposed an Unexposed DPSSC/3 Biological Indicators

Attachment 2

Post-Incubation Pictures






3B Medical

Population Assays on Exposed an Unexposed DPSSC/3 Biological Indicators

Attachment 2 (cont.)





Date: 19 July 2019 Page 8 of 10

Attachment 2 (cont.)





3B Medical

Population Assays on Exposed an Unexposed DPSSC/3 Biological Indicators

Attachment 2 (cont.)

